

Genome assembly strategies

Arturo Vera Ponce de Leon

May 2019





veraponcedeleon.1@osu.edu

History of NGS and Quality control

High throughput sequencing or NGS

High-throughput sequencing

Phase 1: more is better






	2005	GS20	200 000 reads	100 bp
			0.02 Gb/run	
<hr/>				
	2011	GS FLX+	1.2 million reads	750 bp
			0.7 Gb/run	
<hr/>				
	2006	GA	28 million reads	25 bp
			0.7 Gb/run	
<hr/>				
	2011	HiSeq 2000	3 billion reads	2x100 bp
			600 Gb/run	



NORWEGIAN SEQUENCING CENTRE

High-throughput sequencing

Phase 2: smaller is better

1 day		
	→	
0.7 GB/run 700 bp reads		10 hrs GS Junior from Roche/454 0.04 GB/run 400 bp reads
	→	
600 GB/run 2x100 bp reads 10 day		27 hrs MiSeq from Illumina 4.5 GB/run 2x150 bp reads
		
		3 hrs PGM from Ion Torrent/ Life Technologies 0.01, 0.1 or 1 GB/run 100 or 200 bp reads




NORWEGIAN SEQUENCING CENTRE

Sequencing-by-synthesis categories. SBS is a term used to describe numerous DNA-polymerase-dependent methods


454 and IonTorrent sequencing

Template immobilization strategies.




Sequencing by synthesis: single-nucleotide addition approaches.

a 454 pyrosequencing (Roche)





Single nucleotide addition

Only one dNTP species is present during each cycle; multiple identical dNTPs can be incorporated during a cycle, increasing emitted light



b Ion Torrent (Thermo Fisher)




Illumina technology

Template immobilization strategies.

b Solid-phase bridge amplification (Illumina)


Template binding

Free templates hybridize with slide-bound adapters



Bridge amplification

Distal ends of hybridized templates interact with nearby primers where amplification can take place



Cluster generation

After several rounds of amplification, 100–200 million clonal clusters are formed


Patterned flow cell

Microwells on flow cell direct cluster generation, increasing cluster density




Sequencing by synthesis: cyclic reversible termination approaches.

a Illumina




Nucleotide addition

Fluorophore-labelled, terminally blocked nucleotides hybridize to complementary base. Each cluster on a slide can incorporate a different base.



Imaging

Slides are imaged with either two or four laser channels. Each cluster emits a colour corresponding to the base incorporated during this cycle.



Cleavage

Fluorophores are cleaved and washed from flow cells and the 3'-OH group is regenerated. A new cycle begins with the addition of new nucleotides.

High-throughput sequencing

Phase 3: single-molecule



C2 (current) chemistry:
Average read length 2500 bp
36 000 reads
90 MB per 'run'




Real-time sequencing.

(A)




Library preparation comprises the ligation of hairpin adapters (yellow) to double-stranded DNA molecules (blue), thereby creating circular molecules called 'SMRTbells'.

Pacific Biosciences — Real-time sequencing



Nature Reviews | Genetics

(c)



The presence of an epigenetic modification, such as 6-methyladenosine (6 mA), results in a delayed IPD

Trends in Genetics (2018) Vol. 34, No. 9 666-681

Platform Features




Feature	HiSeq2500 - Highoutput	HiSeq2500 – Rapid mode	MiSeq	PacBio RSII
Number of reads	150-180M/lane	100-150M/lane	12-15M (v2) 20-25M (v3)	50-80K/SMRT cell
Read length	2 x 100 bp	2 x 150 bp	2 x 300 bp (v3)	~ 10-20 kb
Yield per lane (PF data)	up to 35 Gb	up to 45Gb	up to 15 Gb	up to 0.4 Gb
Instrument Time	~12-14 days	~2 days	~2 days	~2 hours
Pricing per Gb	\$59 (PE100)	\$53 (PE150)	\$108 (PE300)	\$697

Applications



Platform	454	Illumina HiSeq	Illumina MiSeq*	Ion Torrent	PacBio
resequencing	-	+++	++	-	+
<i>de novo</i>	+++	+	+	+++	+++
metagenomics	+++	++	+	+++	+/-
mRNA	++	+++	++	++	++
miRNA	-	+++	+++	-	-
ChIP	-	+++	++	-	-
DNA meth	-	+++	+	-	!!!
SNP validation	+	-	-	-	++

Multiple technologies diverse features



Which one is the good one?

Yields

A Genome of 1Mb (1 x 10⁶ bases):

- By Sanger:
 - $C = nI/L$
 - $10 = n(500)/1,000,000$
 - $n = 1,000,000 * 10 / 500$
 - 20,000 reads
 - Cost per read ~ 1-2 USD
 - 20,000 USD (~360,000 MX pesos)
- A 454 run ~700Mb (700X)
 - Cost arprox de 20,000 USD
- Un SMRT cell de PacBio (P6-C4) ~150,000 reads (1Gb)
 - Cost 800 USD (~14,400 pesos)
- An Illumina lane ~300 millions of reads (HiSeq2000)
 - An average length of 100 bp = 30 Gb = 30,000 X
 - Cost per lane 2,000 USD (~36,000 pesos)

Coverage:

$$C = nI/L$$

C=Coverage

N=Number of reads

I=Read length

L=Genome size (length) in
bases

30, 000 X coverage ~ \$ 36, 000 Mxpesos

Illumina

10x coverage ~ \$ 360, 000 MX pesos Sanger


**Minimal coverage for SNPs, annotation and
completeness assessments**

> 50 x


Pair end vs Mate Pair

Paired-End Sequencing

(Short-insert paired-end reads)




Mate Pair Sequencing



[Lets watch a very useful video](#)

HTS general analysis flow chart



‘Wet-lab’ experimental design

Bioinformatics hard work

HTS general analysis time flow chart




Fig. 2. Typical whole genome sequencing workflow in a clinical or public health laboratory.

Genome
assembly

Sequence file formats

- Next gen sequence file formats are based on the commonly used

FASTA format

>sequence_ID and optional comments

```
ATTCCGGTGCGGTGCGGTGCTGCCGTGCCGGTGC  
TTCGAAATTGGCGTCAGT
```

- The Phred quality scores per base were added to form the FASTQ format

Sequence file formats

- Illumina Fastq format (fasta format with **Q**uality values for each base)

Read ID

```
@EAS139:136:FC706VJ:2:5:1000:12850 1:Y:18:ATCACG
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA - base calls
+
BBBBCCCC?<A?BC?7@@???????DBBA@@@@A@@ - Base quality+33
```

Space to separate Read

Full read header description

@ <instrument-name>:<run ID>:<flowcell ID>:<lane-number>:<tile-number>: <x-pos>: <y-pos>
<read number>:<is filtered>:<control number>:<barcode sequence>

The phred quality score

Quality score interpretation

$$Q = -10 \log_{10} P \quad \longrightarrow \quad P = 10^{\frac{-Q}{10}}$$

If base quality = 35

$$P = 10^{-35/10} = 0.00032$$

or 1/3200 incorrect

Phred Quality Score	Probability of incorrect base call	Base call accuracy
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1000	99.9%
40	1 in 10000	99.99%
50	1 in 100000	99.999%

```
@EAS139:136:FC706VJ:2:5:1000:12850 1:Y:18:ATCACG
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA - base calls
+
BBBBCCCC?<A?BC?7@@???????DBBA@@@A@@ - Base quality+33
```

ASCII Table

Dec	Hex	Oct	Char	Dec	Hex	Oct	Char	Dec	Hex	Oct	Char	Dec	Hex	Oct	Char
0	0	0		32	20	40	[space]	64	40	100	@	96	60	140	`
1	1	1		33	21	41	!	65	41	101	A	97	61	141	a
2	2	2		34	22	42	"	66	42	102	B	98	62	142	b
3	3	3		35	23	43	#	67	43	103	C	99	63	143	c
4	4	4		36	24	44	\$	68	44	104	D	100	64	144	d
5	5	5		37	25	45	%	69	45	105	E	101	65	145	e
6	6	6		38	26	46	&	70	46	106	F	102	66	146	f
7	7	7		39	27	47	'	71	47	107	G	103	67	147	g
8	8	10		40	28	50	(72	48	110	H	104	68	150	h
9	9	11		41	29	51)	73	49	111	I	105	69	151	i
10	A	12		42	2A	52	*	74	4A	112	J	106	6A	152	j
11	B	13		43	2B	53	+	75	4B	113	K	107	6B	153	k
12	C	14		44	2C	54	,	76	4C	114	L	108	6C	154	l
13	D	15		45	2D	55	-	77	4D	115	M	109	6D	155	m
14	E	16		46	2E	56	.	78	4E	116	N	110	6E	156	n
15	F	17		47	2F	57	/	79	4F	117	O	111	6F	157	o
16	10	20		48	30	60	0	80	50	120	P	112	70	160	p
17	11	21		49	31	61	1	81	51	121	Q	113	71	161	q
18	12	22		50	32	62	2	82	52	122	R	114	72	162	r
19	13	23		51	33	63	3	83	53	123	S	115	73	163	s
20	14	24		52	34	64	4	84	54	124	T	116	74	164	t
21	15	25		53	35	65	5	85	55	125	U	117	75	165	u
22	16	26		54	36	66	6	86	56	126	V	118	76	166	v
23	17	27		55	37	67	7	87	57	127	W	119	77	167	w
24	18	30		56	38	70	8	88	58	130	X	120	78	170	x
25	19	31		57	39	71	9	89	59	131	Y	121	79	171	y
26	1A	32		58	3A	72	:	90	5A	132	Z	122	7A	172	z
27	1B	33		59	3B	73	;	91	5B	133	[123	7B	173	{
28	1C	34		60	3C	74	<	92	5C	134	\	124	7C	174	
29	1D	35		61	3D	75	=	93	5D	135]	125	7D	175	}
30	1E	36		62	3E	76	>	94	5E	136	^	126	7E	176	~
31	1F	37		63	3F	77	?	95	5F	137	_	127	7F	177	

$$\begin{array}{c}
 \text{AAAAA} \\
 + \\
 \text{BBBBB} \\
 \downarrow \downarrow \downarrow \downarrow \downarrow \\
 \text{ASCII val} \quad 66 \quad 67 \\
 \quad \quad \quad -33 \quad -33 \\
 \\
 \text{Q value} \quad \frac{33}{34}
 \end{array}$$


$$Q = -10 \log_{10} P \quad \longrightarrow \quad P = 10^{\frac{-Q}{10}}$$

$$\begin{aligned}
 &> 10^{(-33/10)} \\
 &[1] \ 0.0005011872 = 1/5000 \\
 &> 10^{(-34/10)} \\
 &[1] \ 0.0003981072 = 1/39000
 \end{aligned}$$

Let's play with FasQC to quality control visualization


Open FastQC program

Open in browser:
fastqc_report.html



Quality filter trim galore

Before trimGalore



After trimGalore

